LIPOGEMS PROCESS DEVICE

LGD 60

REF (LGD 60 O, LGD 60 B, LGD 60 S, LGD 60 M)

LGD 120

REF (LGD 120 O, LGD 120 B, LGD 120 S, LGD 120 M)

LGD 240

REF (LGD 240 O, LGD 240 B, LGD 240 S, LGD 240 M)

EN Lipoaspirate fat processing device	PL Urządzenie do przetwarzania tkanki tłuszczowej pobranej metodą lipoaspiracji
IT Dispositivo per la processazione del tessuto adiposo lipoaspirato	ZH (简体中文) 脂肪抽吸物处理器械
DE Aufbereitungsgerät für abgesaugtes Fett	AR (العربية) جهاز معالجة شفط الدهون
ES Producto para procesado de grasa lipoaspirada	JA 吸引脂肪組織処理装置
FR Dispositif de traitement des graisses par lipoaspiration	TR Lipoaspirat yağ işleme cihazı
PT Dispositivo de processamento de gordura lipoaspirada	DA Behandlingsudstyr til fedtsugning
HU Lipoaspirátum (leszívott zsírszövet) feldolgozó eszköz	SV Lipoaspiratprocessenhet
HR Uređaj za obradu masti lipoaspirata	RO Dispozitiv de procesare a grăsimilor lipoaspirate
EL Συσκευή επεξεργασίας λίπους λιποαναρρόφησης	LT Lipoaspirato riebalų apdorojimo įrenginys

Indications for Use / Intended Use

The Lipogems System is a sterile medical device intended for the closed-loop processing of lipoaspirate tissue in medical procedures involving the harvesting, concentrating and transferring of autologous adipose tissue harvested with a legally marketed lipoplasty system. The device is intended for use in the following surgical specialties when the transfer of harvested adipose tissue is desired: orthopedic surgery, arthroscopic surgery, neurosurgery, gastrointestinal and affiliated organ surgery, urological surgery, general surgery, gynecological surgery, thoracic surgery, laparoscopic surgery, and plastic and reconstructive surgery when aesthetic body contouring is desired. Only legally marketed accessory items, such as syringes, should be used with the system. If harvested fat is to be transferred, the harvested fat is only to be used without any additional manipulation.

Device description

The device comes pre-assembled and consists of:

- a processing unit, which contains stainless steel balls and sieves,
- a saline input line, connected to the blue/orange/yellow/pink side (Figure 1, ref. 7) of the processing unit,
- an access port with Luer-Lock connection and self occluding valve (Figure 1, ref. 6) at blue/orange/yellow/pink side, to load material to be processed,
- a drainage line, connected to the gray side of the processing unit,
- an access port with Luer-Lock connection and self-occluding valve (Figure 1, ref.10) at gray side, to unload the processed lipoaspirate, and
- a collection bag for waste fluid (Figure 1, ref. 16).

None of the Lipogems System components are manufactured using natural rubber latex or phthalates.

Principle of operation

The processing unit, filled with normal saline solution is maintained in a continuous flow condition by gravity. After saline solution priming, adipose tissue inserted in the device undergoes a first treatment to reduce and micronize lipid clusters by means of an input sieve. The mechanical action exerted by the balls, obtained by shaking the processing unit, allows emulsion of the lipid mass and consequent reduction of the clusters. This processing diminishes tissue disaggregation, thus preserving the cell and tissue micro-architecture of the adipose tissue. The continuous flow of normal saline solution eliminates residues of oil emulsion and any remaining blood components.

Packaging and Sterility

The device is sterilized by means of ethylene oxide. Before use, carefully check the integrity of packaging to make sure that it has not suffered any damage, which may

compromise the sterility of its contents. After opening the package, it is absolutely necessary, and the operator's responsibility, to apply a strictly aseptic technique. Any error in the handling or during transfer to the sterile field may jeopardize the sterility of the device, of the surgical intervention, and pose a serious risk of complications for the patient, such as infection and sepsis.

Instructions for use

1. Open the primary package and remove the device.
2. Verify the connection between the main unit and the drain collection bag (Figure 1, ref. 14 and 16) and place the drain collection bag on or near the ground. Prior to using the device, make sure that the valves (Figure 1, ref. 6 and 10) are secured firmly to the processing unit's gray and colored end caps, exercising care not to over tighten, or connection breakage may occur.
3. Connect the terminal spike of the input line (Figure 1, ref. 1) to the saline bag by piercing the bag connection port (recommended: 2000 ml bag for LGD 60, 3000 ml for LGD 120 and 5000 ml for LGD 240). Then hang the bag at a higher position than the processing unit (Figure 1, ref. 8). The spike is equipped with a soft dripping chamber (Figure 1, ref. 2), which allows visualization of saline flow and, eventually, to carry out replacement of the saline bag without introducing air into the circuit.
4. Verify that all five clamps (Figure 1, ref. 3, 5, 11, 13 and 15) in the circuit are open. Place the device in vertical position with the gray side (Figure 1, ref. 9) upward and the blue/orange/yellow/pink side downward. The saline solution fills the unit and flows into the collection bag (Figure 1, ref. 16). Check if there are air bubbles in the processing unit and, if so, remove them manually by orienting and shaking the processing unit until complete air removal. Verify that there are no impediments to the flow of saline.
5. Prior to loading the device with lipoaspirate, make sure that the end cap valves (Figure 1, ref. 6 and 10) are secured firmly to the process unit caps. Do not over-tighten, otherwise connection breakage may occur. Keep the processing unit (Figure 1, ref. 8) in a vertical position at eye level, with the blue/orange/yellow/pink end cap (Figure 1, ref. 7) upwards. Close the clamp very near to the blue/orange/yellow/pink head input line (Figure 1, ref. 5), and prepare to fill the device with lipoaspirate by connecting the syringe to the self-occluding valve (Figure 1, ref. 6) on the blue/orange/yellow/pink end cap, while ensuring that the syringe is directly in line with vertical axis of the valve, and not applying any lateral or off-axis force. Insert the lipoaspirate in the processing unit (Figure 1, ref. 8) by depressing the syringe plunger. Repeat this step until the desired amount of lipoaspirate is injected, introducing the tissue within the limits of 15 cc at a time for LGD60, 40 cc at a time for LGD 120 and 80 cc at a time for

LGD 240. Do not exceed these indicated quantities of lipoaspirate.

6. Open the blue/orange/yellow/pink end cap clamp (Figure 1, ref. 5) and restore saline flow until the solution exiting the processing unit appears clear.
7. Before further processing the adipose tissue, ensure that the input clamp and drain line clamp (Figure 1, ref. 5 and 11) are closed very near to the blue/orange/yellow/pink and gray end caps, respectively. Process the adipose tissue by vigorously shaking the processing unit (Figure 1, ref. 8) for 30 seconds. Next, open the blue/orange/yellow/pink and gray end cap clamps (Figure 1, ref. 5 and ref. 11) and restore saline flow for approximately 10 seconds. Repeat this cycle of shaking and saline flow for at least 2 minutes, until the adipose tissue appears compact, light yellow in color, and floats on the transparent saline solution.
8. Position the processing unit (Figure 1, ref. 8) vertically with the gray end cap upwards (Figure 1, ref. 9). The processed tissue will float to the top of the processing unit. Close the clamp located on the drain line close to the gray end cap (Figure 1, ref. 11). Connect a 10 ml luer lock syringe at each end cap valve of the processing unit. (Figure 1, ref. 6 and 10). Before connecting, make sure that the syringe plungers are completely inserted, leaving no air in the syringe.
9. Withdraw approximately 10 ml of saline from the input line using the syringe connected at the blue/orange/yellow/pink side (Figure 1, ref. 6). If during this phase some fibrous tissue enters the syringe, detach the syringe, discard the fibrous tissue and fill the syringe again. Next, close the clamp (Figure 1, ref. 5) located on the input line. With both end clamps closed (Figure 1, ref. 5 and 11), firmly reinject the saline into the processing unit. The syringe connected to the valve on the grey end cap (Figure 1, ref. 10), now on the upward side of the processing unit, will be consequently filled with the processed tissue as the saline is injected at the opposite end of the processing unit. Do not attempt to re-inject saline solution with a syringe other than the recommended 10 ml luer lock. Firmly depress the plunger, but avoid over-forcing. Do not use other clamps, clamping devices or other instruments to occlude input and output lines.
10. Remove the syringe with the processed adipose tissue. If desired, the user may obtain more tissue by repeating step 9.
11. If still additional tissue is desired, obtain a new lipoaspirate, and then process and collect by repeating steps 5 to 10. Only one repeat loading and processing cycle is permitted, as the device may not function effectively for subsequent processing.
12. When finished, close the clamps near the saline bag and collection bag (Figure 1, ref. 3 and 13), disconnect the processing device and dispose of all components in accordance with local protocols.

Storage

The device must be kept in its original unopened packaging, in a dry, clean, indoor environment, at a temperature not exceeding 50°C, away from heat sources and avoiding direct and prolonged exposure to sunlight.

Duration, sterility, disposal

The expiry date is stated on the label of the primary package. Do not use the device after the expiry date. The device must not be used if the primary packaging is not present and fully intact.

The unused content of open packages must not be used and must be disposed of as waste. LGD 60, LGD 120 and LGD 240 are sterilized by ethylene oxide and must not be re-sterilized. Dispose of used, expired, unusable or damaged devices by following local protocols applicable to the specific type of hospital waste.

Contraindications

The device (LGD 60, LGD 120, LGD 240) must be used only by licensed and trained physicians who are aware of the potential side effects, complications, limitations and contraindications to harvesting procedures of subcutaneous adipose tissue, tissue processing and subsequent grafting. It is not recommended to perform the procedure on patients suffering from coagulation disorders or in patients under anticoagulant therapy.

Warnings

- This device will not, in and of itself, produce significant weight reduction.
- This device should be used with extreme caution in patients with chronic medical conditions, such as diabetes; heart, lung, or circulatory system disease; or obesity.
- The volume of blood loss and endogenous body fluid loss may adversely affect intra and/or postoperative hemodynamic stability and patient safety. The capability of providing adequate, timely replacement is essential for patient safety.

Precautions

- This device is designed to remove localized deposits of excess fat through small incisions and subsequently transfer the tissue back to the patient.
- Results of this procedure will vary depending upon patient age, surgical site, and experience of the physician.
- The amount of adipose tissue removed should be limited to that necessary to achieve the desired cosmetic effect.
- Results of this procedure may or may not be permanent.
- The use of the device is restricted to licensed physicians, trained to perform procedures for harvesting, processing and grafting subcutaneous adipose tissue.

- Single use device. Do not reuse. The device cannot and must not be re-sterilized or reused. This device has been designed and produced for a single use. Reuse may result in contamination, loss of mechanical integrity and malfunction leading to loss of mechanical functionality and harm to the patient.
- When loading, do not exceed the maximum volume of tissue, which is 15 ml at a time for LGD 60, 40 ml at a time for LGD 120 and 80 ml at a time for LGD 240 for a maximum of 2 times in order to avoid sieve obstruction or other malfunctioning.
- Ensure all components of the device are present. The manufacturer assumes no responsibility for damages to persons or property if non original components are used.

- After use, the device poses a contamination risk, even if cleansed or disinfected. Avoid accidental direct or indirect contact with contaminated devices.
- Do not use the LGD 60, the LGD 120 and LGD 240 for any other purpose than that detailed in the Indications for Use/Intended Use statement.

The user is required to report any serious incident that has occurred, in relation to the device, to the manufacturer and to the Competent Authority of the Member State in which the user is established in the manner prescribed by this Authority and by the Medical Devices Regulation.

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Lipogems International SpA